

Certificate of Compliance



We confirm that all parts in contact with the dissolution media are made of inert material such:

2830 / 9353 Basket without Lock	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462
2985 Paddle with Suction Opening	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462
5301 Basket Holder with Suction Opening	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462
7525 Paddle with Stirring Shaft	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462
7558 Basket Holder with Stirring Shaft	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462
7710 Suction Head with Filter Insert	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462

Aesch, 2019-10-14

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland