



Part Paddle with Suction Opening
 Part Number 02985-2
 Serial Number N13003794 – N13003941



Test Certificate

				USP / Ph.Eur. Tolerances
Dimensions	Paddle	Height	19.0 mm	± 0.5 mm
		Upper edge length	74.5 mm	± 0.5 mm
		Lower edge length	42.0 mm	± 1.0 mm
		Thickness	4.0 mm	± 1.0 mm
		Radius	1.2 mm	± 0.2 mm
		Radius	41.5 mm	± 1.0 mm
	Shaft	Diameter	9.4 – 10.1 mm	--
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2024-02-24

Tayyib Özüdogru
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland