



Part Basket Holder with Stirring Shaft
Part Number 07558-07
Serial Number P13 000501 – P13 001084



			USP / Ph.Eur. Tolerances
Dimensions	Shaft Ø	9.4 – 10.1 mm	--
	Cap Ø	25.0 mm	± 3.0 mm
	Cap size	5.1 mm	± 0.5 mm
	Clear opening	20.2 mm	± 1.0 mm
	Vent hole Ø	2.0 mm	± 0.5 mm
Material	Stainless Steel (electrolytic polished)		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2022-10-11

mp

11-10-2022 09:31:12

Moreshwar Pulekar
SOTAX Quality Management
C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India

Test Certificate