



Test Certificate

Part Suction Head with filter insert
Part Number 07710-1
Serial Number N15004000 - N15004750



				USP / Ph.Eur. Tolerances
Dimensions	Shaft	Diameter	10 mm	9.4 – 10.1mm
Material	Suction head PVC			
	Shaft Stainless Steel (electrolytic polished)			
	Suction tube PTFE			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2024-08-23

Michael Schmidlin
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland