



Part

Paddle with Stirring Shaft

Part Number

08365

Serial Number

N23000184 - N23000195



Test Certificate

| | | | | USP / Ph.Eur. Tolerances |
|------------|---|-------------------|---------------|-----------------------------|
| Dimensions | Paddle | Height | 19.0 mm | ± 0.5 mm |
| | | Upper edge length | 74.5 mm | ± 0.5 mm |
| | | Lower edge length | 42.0 mm | ± 1.0 mm |
| | | Thickness | 4.0 mm | ± 1.0 mm |
| | | Radius | 1.2 mm | ± 0.2 mm |
| | | Radius | 41.5 mm | ± 1.0 mm |
| | Shaft | Diameter | 9.4 – 10.1 mm | -- |
| Material | Stainless Steel (electrolytic polished) | | | |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2025-04-10

Buner Cornelius
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland