

Test Certificate

Part: **Basket, 40-Mesh SS, Dissolution**
 Number: 63127
 Serial Number: QLA86449-QLA86590

For SOTAX AT - Dissolution Testing Systems



		USP Tolerances	
Dimensions:	Basket rim O.D. Ø	25.0 mm	± 3.0 mm
	Height	37.0 mm	± 3.0 mm
	Wire Ø	0.25 mm	0.22 – 0.31 mm
	Wire openings	0.4 mm	0.36 - 0.44 mm
	Screen opening height	27.0 mm	± 1.0 mm
	Screen O.D. Ø	22.2 mm	± 1.0 mm
	Clear opening Ø	20.2 mm	± 1.0 mm

Material: **Stainless Steel 316**

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP Pharmacopeia requirements (at time of delivery).

A handwritten signature in blue ink that reads "Ajay".

SOTAX Quality Management
 Ajay Pydah
 July 23, 2018